

Comprehensive Histopathology

<https://dmm.biologists.org/content/7/5/515.long>

This was a collaboration with the Centre for Modeling Human Disease, Toronto Centre for Phenogenomics, Canada:

For select lines, the complete tissue blocks from our biobank were sent off. These were sectioned at 4 μ m for routine hematoxylin and eosin (H&E) staining prior to being screened by veterinary pathologists. Images were captured using a microscope-mounted Olympus DP71 digital camera (Olympus Life Science Imaging Systems Inc., Markham, ON, Canada).

The age at which the mice were culled and the diet they were fed was pipeline dependent.