



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	0.6232	1.462E-03
24S/3S	-0.9297	3.724E-06
A/24S	-0.1309	5.620E-01
24S/C	-0.3065	6.930E-02

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
58.6265	53.9702	113.7631	50.2462