



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	0.1058	6.191E-01
24S/3S	-0.6521	1.890E-03
A/24S	-0.3402	1.229E-01
24S/C	-0.5463	9.362E-04

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
46.6673	52.6356	68.9819	29.32