



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	-1.3113	4.912E-02
24S/3S	-1.8793	1.772E-02
A/24S	2.7556	3.459E-04
24S/C	-3.1946	5.860E-08

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
1.2267385	0.299293	0.464121	0.983362