



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	1.9420	3.681E-03
24S/3S	-0.3531	7.070E-01
A/24S	-2.9531	6.807E-03
24S/C	1.5925	6.555E-03

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
0.373729	0.638253	1.90574	0.179235