



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	-0.8767	5.892E-04
24S/3S	-0.3130	2.788E-01
A/24S	1.1841	3.505E-07
24S/C	-1.1897	1.551E-09

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
195.1565	317.474	151.8075	225.69