



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	1.9834	1.642E-04
24S/3S	-0.3150	6.168E-01
A/24S	-1.3679	4.617E-02
24S/C	1.6706	2.982E-04

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
0.78022	0.827085	4.77365	1.10712