



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	-0.6276	1.405E-02
24S/3S	-0.0497	8.899E-01
A/24S	0.6678	6.633E-03
24S/C	-0.6774	5.978E-04

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
69.10675	83.3028	72.3192	78.8861