



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	-0.6693	1.790E-03
24S/3S	0.1352	6.172E-01
A/24S	0.6358	1.273E-03
24S/C	-0.5341	1.406E-03

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
87.0106	66.7372	99.3138	109.116