



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	0.8437	2.792E-04
24S/3S	-0.2308	4.123E-01
A/24S	-2.7817	3.157E-22
24S/C	0.6129	1.634E-03

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
79.8773	158.886	234.399	20.819