



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	-0.4133	4.732E-02
24S/3S	-0.3184	1.695E-01
A/24S	0.5054	1.223E-02
24S/C	-0.7317	5.372E-06

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
72.62465	83.487	78.7167	71.4382