



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	0.1817	4.850E-01
24S/3S	0.5528	1.012E-01
A/24S	-0.4612	1.398E-01
24S/C	0.7345	1.010E-03

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
19.1942	32.779	52.7246	27.4404