



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	0.6407	1.280E-03
24S/3S	0.1484	5.688E-01
A/24S	-0.0662	7.812E-01
24S/C	0.7891	1.129E-06

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
37.8475	63.6805	105.42065	71.9949