



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	1.9596	1.449E-03
24S/3S	-1.1428	6.548E-02
A/24S	-3.3481	1.214E-03
24S/C	0.8183	1.835E-01

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
2.1469	6.21936	9.382645	0.384829