



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	-0.6070	1.394E-02
24S/3S	-0.8096	1.427E-03
A/24S	0.8115	6.975E-04
24S/C	-1.4167	2.274E-13

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
150.5985	198.672	121.90395	112.807