



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	0.4639	5.660E-02
24S/3S	-0.6197	1.856E-02
A/24S	-0.6082	1.755E-02
24S/C	-0.1558	5.118E-01

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
51.90055	48.1005	94.324	35.5346