



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	-0.1542	4.948E-01
24S/3S	-0.6045	3.752E-03
A/24S	-0.6256	2.727E-03
24S/C	-0.7587	1.965E-06

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
56.81335	35.2935	71.2276	25.1625