



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	5.0166	2.383E-62
24S/3S	-1.2826	6.167E-06
A/24S	0.6498	2.552E-02
24S/C	3.7355	3.135E-42

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
0.807991	4.84726	29.4646	19.7871