



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	2.2933	5.437E-07
24S/3S	-0.7386	1.492E-01
A/24S	-2.4429	2.231E-04
24S/C	1.5561	1.017E-04

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
1.283995	0.967121	8.61816	0.830741