



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	-0.7828	1.957E-04
24S/3S	0.3615	1.199E-01
A/24S	0.5733	3.121E-03
24S/C	-0.4213	1.162E-02

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
428.198	422.851	491.631	555.105