



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	-0.7519	3.675E-04
24S/3S	0.3716	1.138E-01
A/24S	-0.5199	1.431E-02
24S/C	-0.3803	2.298E-02

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
31.2904	32.1327	37.0292	19.5032