



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	3.2484	4.465E-07
24S/3S	-1.8755	4.443E-03
A/24S	-3.7505	4.196E-03
24S/C	1.3768	1.829E-02

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
0.42503385	2.49084	5.084345	0.101036