



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	-1.3854	3.850E-11
24S/3S	0.4922	3.397E-02
A/24S	0.1534	4.643E-01
24S/C	-0.8932	1.713E-08

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
91.7115	45.4597	73.9381	62.9803