



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	0.9856	3.197E-07
24S/3S	-0.6515	1.796E-03
A/24S	-1.2679	1.970E-09
24S/C	0.3341	4.466E-02

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
59.8744	91.9441	159.7785	36.1959