



recommended PAdjust cut-off < 0.01

| <b>comp</b> | <b>log2FC</b> | <b>PAdjust</b> |
|-------------|---------------|----------------|
| 3S/C        | -1.7612       | 7.139E-07      |
| 24S/3S      | 0.8923        | 1.581E-02      |
| A/24S       | 1.1274        | 4.012E-04      |
| 24S/C       | -0.8691       | 1.448E-03      |

RPKM values are shown below - should not be compared among samples!

| <b>cerc</b> | <b>3Hsomule</b> | <b>24Hsomule</b> | <b>adult</b> |
|-------------|-----------------|------------------|--------------|
| 1.69337     | 1.1943          | 1.1054035        | 2.38123      |