



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	0.8668	1.926E-03
24S/3S	0.4939	1.280E-01
A/24S	-0.9151	2.142E-03
24S/C	1.3608	2.217E-09

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
16.189	25.3505	56.6436	25.4503