



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	0.1621	4.434E-01
24S/3S	0.1481	5.761E-01
A/24S	-0.8426	7.085E-05
24S/C	0.3101	5.942E-02

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
30.16545	29.8796	61.15065	24.2402