



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	1.7472	3.659E-09
24S/3S	-0.9034	5.967E-03
A/24S	-2.3718	8.545E-10
24S/C	0.8438	8.913E-04

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
19.7273	28.8687	85.1637	8.06559