



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	1.3034	2.377E-06
24S/3S	-0.6163	4.010E-02
A/24S	-0.7653	1.212E-02
24S/C	0.6871	4.801E-03

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
12.47725	35.3191	42.1378	13.3696