



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	-0.6231	3.300E-03
24S/3S	0.6126	5.384E-03
A/24S	-1.5189	5.880E-13
24S/C	-0.0105	9.804E-01

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
32.1682	23.6334	43.6812	12.9857