



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	-0.7636	3.099E-04
24S/3S	0.4278	6.164E-02
A/24S	0.0861	6.977E-01
24S/C	-0.3358	5.190E-02

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
170.0145	150.169	199.059	166.885