



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	0.6368	1.190E-01
24S/3S	0.1443	8.062E-01
A/24S	-1.2841	5.692E-03
24S/C	0.7813	1.689E-02

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
0.905156	1.26971	2.50697	0.737433