



recommended PAdjust cut-off < 0.01

<b>comp</b>	<b>log2FC</b>	<b>PAdjust</b>
3S/C	-0.1256	6.140E-01
24S/3S	-0.5299	1.408E-02
A/24S	-0.1356	5.603E-01
24S/C	-0.6556	6.958E-05

RPKM values are shown below - should not be compared among samples!

<b>cerc</b>	<b>3Hsomule</b>	<b>24Hsomule</b>	<b>adult</b>
97.6666	149.261	126.876	65.6222