



recommended PAdjust cut-off < 0.01

| <b>comp</b> | <b>log2FC</b> | <b>PAdjust</b> |
|-------------|---------------|----------------|
| 3S/C        | -0.1154       | 8.525E-01      |
| 24S/3S      | -0.7790       | 6.673E-02      |
| A/24S       | -1.9852       | 3.965E-04      |
| 24S/C       | -0.8945       | 6.380E-03      |

RPKM values are shown below - should not be compared among samples!

| <b>cerc</b> | <b>3Hsomule</b> | <b>24Hsomule</b> | <b>adult</b> |
|-------------|-----------------|------------------|--------------|
| 6.06446     | 5.69939         | 7.406545         | 0.960307     |